

Product datasheet for **MG227530**

Trp63 (NM_001127259) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trp63 (NM_001127259) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Trp63
Synonyms:	AI462811; delta; Ket; p6; p7; p51/p; P51/P63; P63; P73I; TAp; Tp63; Trp5; Trp53rp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG227530 representing NM_001127259
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATTTTGAACCTTCACGGTGTGCCACCTACAGTACTGCCCGACCCTTACATCCAGCGTTTCATAG
 AAACCCAGCTCATTCTCGTGGAAAGAAAGTTATTACAGATCTGCCATGTCGCAGAGCACCCAGACAAG
 CGAGTTCTCAGCCAGAGGCTTCCAGCATATCTGGATTTTCTGGAACAGCCTATATGCTCAGTACAG
 CCCATCGAGTTGAACTTTGTGGATGAACCTTCCGAAAATGGTGCAACAAACAAGATTGAGATTAGCATGG
 ATTGTATCCGCATGCAAGACTCAGACCTCAGTGACCCCATGTGGCCACAGTACACGAACCTGGGGCTCCT
 GAACAGCATGGACCAGCAGATTGAGAACGGCTCCTCGTCCACCAGCCCTACAACACAGACCACGCACAG
 AATAGCGTGACGGCGCCCTCGCCCTATGCACAGCCAGCTCCACCTTTGATGCCCTCTCTCCATCCCCTG
 CCATTCCCTCAAACACAGATTACCCGGGCCACACAGCTTCGATGTGCCTTCCAGCAGTCAAGCACTGC
 CAAGTCAGCCACCTGGACGTATCCACCGAACTGAAGAAGCTGTACTGCCAGATTGCGAAGACATGCCCC
 ATCCAGATCAAGGTGATGACCCACCCACAGGGCGCTGTTATCCGTGCCATGCCTGTCTACAAGAAAG
 CTGAGCATGTACCCGAGGTTGTGAAACGATGCCCTAACCATGAGCTGAGCCGTGAGTTCAATGAGGGACA
 GATTGCCCTCCCAGTCACTGATTGAGTAGAAGGGAACAGCCATGCCAGTATGTAGAAGATCCTATC
 ACGGGAAGGCAGAGCGTGCTGGTCCCTTATGAGCCACCACAGGTTGGCACTGAATTCACAACAGTCTGT
 ACAATTTTCATGTGAACAGCAGCTGCGTCGGAGGAATGAACAGACGTCCAATTTTAAATCATCGTTACTCT
 GGAAACCAGAGATGGGCAAGTCTGGGCCGACGGTGCTTTGAGGCCCGGATCTGTGCTTGCCAGGAAGA
 GACCGGAAGGCAGATGAAGACAGCATCAGAAAGCAGCAAGTATCGGACAGCGCAAAGAACGGCGATGGTA
 CGAAGCGCCCTTCCGTGAGAATACACACGGAATCCAGATGACTTCCATCAAGAAACGGAGATCCCAGA
 TGATGAGCTGCTGTACCTACCAGTGAGAGGTCGTGAGACGTACGAGATGTTGCTGAAGATCAAAGAGTCA
 CTGGAGCTCATGAGTACCTCCCTCAGCACAGTACGAAACGTACAGGCAGCAGCAGCAGCAGCAGCACC
 AGCACCTACTTCAGAAACAGACCTCGATGCAGTCTCAGTCTTATATGGCAACAGTTCACCCACCTGAA
 CAAAATGAACAGCATGAACAAGTGCCTTCCGTGAGCCAGCTTATCAACCCACAGCAGCGCAATGCCCTC
 ACTCCCACCACCATGCCTGAGGGCATGGGAGCCAAACATTCCTATGATGGCACTCACATGCCAATGGCTG
 GAGACATGAATGGACTCAGCCCTACCAAGCTCTCCCTCCTCCACTCTCCATGCCCTCCACCTCCCCTG
 CACCCACACCCGCTACCCACAGACTGCAGCATTGTCAGTTTCTTAGCAAGGTTGGGCTGCTCATCA
 TGCTTGACTATTTACGACCCAGGGGCTGACCACCATCTATCAGATTGAGCATTACTCCATGGATGATT
 TGGCAAGTCTGAAGATCCCTGAACAGTTCGACATGCCATCTGGAAGGGCATCCTGGACCACAGGCAGCT
 GCACGACTTCTCCTCACCTCCTCATCTCCTGAGGACCCCAAGTGGTGCCTCTACCGTCAGTGTGGGCTCC
 AGTGAGACCCGTGGTGAACGTGTGATCGATGCCGTGCGCTTTACCTCCGCCAGACCATCTCTTTCCAC
 CCCGTGACGAGTGGAAATGATTTCAACTTTGACATGGATTCTCGTCGCAACAAGCAGCAGCGTATCAAAGA
 GGAAGGAGAA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG227530 representing NM_001127259
 Red=Cloning site Green=Tags(s)

MNFETSRCATLQYCPDPYIQRFIETPAHFSWKESYYRSAMSQSTQTSEFLSPEVFQHIWDFLEQPICSVQ
 PIELNFVDEPSENGATNKIEISMDCIRMQDSDLSDPMWPQYTNLGLLNSMDQQIQNGSSSTSPYNTDHAQ
 NSVTAPSPYAQPSSTFDALSPSPAIPSNTDYPGPHSFDVSFQSSSTAKSATWTYSTELKKLYCQIAKTCP
 IQIKVMTPPPQGAVIRAMPVYKKAIEHVTEVVKRCPNHEL SREFNEGQIAPPSHLIRVEGNSHAQYVEDPI
 TGRQSVLVPYEPPQVGTEFTTVLYNFMCNSSCVGGMNRRPILIIIVTLETRDQVLRRCFEARICACPGR
 DRKADEDSIRKQVSDSAKNGDGTKRPFQNTHTGIQMTSICKRRSPDDELLYLPVRGRETYEMLLIKES
 LELMQYLPQHTIETRQQQQQHQLLQKQTSMQSSSYGNSSPPLNKMNSMKNLPSVSQLINPQQRNAL
 TPTTMEPGMANIPMMGTHMPMAGDMNGLSPTQALPPPLSMPSTSHCTPPPPYPTDCSIVSFLARLGCSS
 CLDYFTTQGLTTIYQIEHYSMDDLASLKIPEQFRHAIWKILDHRQLHDFSSPHLLRTPSGASTVSVGS
 SETRGERVIDAVRFTLRQTISFPPRDEWDFNFDMDSRRNKQQRKEEGE



TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja1796_c04.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_001127259

ORF Size: 2040 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127259.1](#), [NP_001120731.1](#)

RefSeq Size: 5382 bp

RefSeq ORF: 2043 bp


Locus ID: 22061

UniProt ID: [O88898](#)

Cytogenetics: 16 17.37 cM

Gene Summary:

This gene encodes tumor protein p63, a member of the p53 family of transcription factors involved in cellular responses to stress and development. The family members include tumor proteins p53, p63, and p73, which have high sequence similarity to one another. This similarity allows p63 and p73 to transactivate p53-responsive genes causing cell cycle arrest and apoptosis. The family members can interact with each other in many ways, including direct and indirect protein interactions. This results in mutual regulation of target gene promoters. Tumor protein p63 $-/-$ mice have several developmental defects which include the lack of limbs and other tissues, such as teeth and mammary glands, which develop as a result of interactions between mesenchyme and epithelium. Both alternative splicing and the use of alternative promoters result in multiple transcript variants encoding different protein isoforms.[provided by RefSeq, Dec 2009]

Product images:

Circular map for MG227530