

Product datasheet for **MG227543**

Dcx (NM_001110223) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dcx (NM_001110223) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Dcx
Synonyms: Dbct
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG227543 representing NM_001110223
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAACCTGATTTTGGACATTTTGACGAACGAGACAAAGCATCTAGAAATATGAGAGGGTCACGGATGA
ATGGACTCCAAGTCCCCTCATAGTGGCCACTGTAGCTTCTACAGAACCAGAACCTTGCAGGCATTAAG
TAATGAGAAGAAGGCCAAGAAGGTACGTTTCTACCGCAATGGGGACCGTTACTTCAAGGGGATTGTGTAC
GCTGTTTCTTCTGACCGTTTTCTGAGTTTTGATGCGTTGCTGGCTGACCTGACCCGATCCTTGTCTGACA
ACATTAACCTGCCTCAGGGAGTGCCTACATTTATACCATTGACGGATCCAGGAAGATTGGAAGCATGGA
TGAAGTGAAGAAGGGGAAAGCTATGTCTGCTCCTCAGACAACTCTTTAAAAAGTTGAGTACACCAAG
AATGTCAACCCCAACTGGTCTGTCAACGTAAGACATCTGCCAACATGAAAGCCCCCAGTCCCTTGCTA
GCAGCAACAGTGTCAAGCCAGAGAGAACAAGGACTTTGTGCGCCCCAACTTGTGACCATCATTGCGAG
CGGGGTGAAGCCACGGAAGGCTGTGCGCGTCTTCTCAACAAGAAAACAGCCCACTTTTCGAGCAGGTC
CTGACTGACATCACAGAAGCGATCAAACCTGAAACCGGAGTTGTCAAAAACTCTACACCCTTGATGGAA
AGCAGGTCACCTGTCTCCATGATTTCTTGGTGATGATGATGTGTTTATTGCTTGTGGTCTGAAAAATT
CCGCTATGCTCAAGATGATTTCTCCTTGGATGAGAATGAATGCAGAGTCATGAAAGGGAATCCATCTGCC
GCAGCTGGCCAAAGGCTTCCCCAACCTCAAAGACATCTGCTAAAAGCCCAGGCCAATGCGCCGCA
GCAAGTCTCCAGCTGACTCAGGTAACGACCAAGACGCAAATGAACTTCCAGCAGTCAGCTCAACACC
TAAGTCAAAGCAGTCTCCTATCTCTACACCCACAAGCCCTGGAAGTCTGCGGAAGCACAAGGACCTGTAC
CTGCCGCTGTCATTGGATGACTCTGATTCACTTGGCGATTCCATG

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG227543 representing NM_001110223
Red=Cloning site Green=Tags(s)

MELDFGHFDERDKASRNMRGSRMNGLPSPTHSAHCSFYRTRTLQALSNEKKAKKRVFYRNGDRYFKGIVY
 AVSSDRFRSFDALLADLTRSLSDNINLPQGVRYIYIDGSRKIGSMDELEEGESYVCSSDNFFKKVEYTK
 NVNPNWSVNVKTSANMKAPQSLASSNSAQARENKDFVRPKLVTIIRSGVKPRKAVRVLLNKKTAHSFEQV
 LTDITEAIKLETGVVKKLYTLDGKQVTLHDFFGDDDFIACGPEKFRYAQDDFSLDENECRVMKGNPSA
 AAGPKASPTPKTSAKSPGPMRRSKSPADSGNDQDANGTSSSQLSTPKSKQSPISPTSPGSLRKHKDL
 LPLSLDSDSLGDSM

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001110223

ORF Size: 1098 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001110223.1](#), [NP_001103693.1](#)

RefSeq Size: 9000 bp

RefSeq ORF: 1101 bp


Locus ID: 13193

UniProt ID: [O88809](#)

Cytogenetics: X F2

Gene Summary: This gene encodes a member of the doublecortin family. The protein encoded by this gene is a cytoplasmic protein and contains two doublecortin domains, which bind microtubules. In the developing cortex, cortical neurons must migrate over long distances to reach the site of their final differentiation. The encoded protein appears to direct neuronal migration by regulating the organization and stability of microtubules. In addition, the encoded protein interacts with LIS1, the regulatory gamma subunit of platelet activating factor acetylhydrolase. Studies in knockout mice lacking this gene and the LIS1 gene suggest that the molecular interaction of these two genes is important in both in neuronal migration and neurogenesis, and there is a cortical role of this gene in nuclear translocation and positioning of the mitotic spindle in radial glial mitotic division. Multiple transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Sep 2010]

Product images:



Circular map for MG227543