

Product datasheet for **MK200005**

Nr4a1 Mouse qPCR Template Standard (NM_010444)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_010444

Synonyms: Gfrp; GFRP1; Hbr-1; Hbr1; Hmr; N10; NGFI-B; NGFIB; NP10; nur77; TIS1; TR3

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP209119) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)