

Product datasheet for **MK200008**

Unc119 Mouse qPCR Template Standard (NM_011676)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_011676
Synonyms:	D11Bhm52; HRG; HRG4; MRG; MRG4; Rg; Rg4; Rtg; Rtg4; UNC11; Unc119h
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP220997) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

