

Product datasheet for **MK200088**

Pdia3 Mouse qPCR Template Standard (NM_007952)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007952
Synonyms:	58kDa; Erp; ERp57; ERp60; ERp61; Grp58; PDI; PDI-Q2; PI-PLC; Plca; PLC[a]
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP211917) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)