

Product datasheet for **MK200113**

Cdkn1a Mouse qPCR Template Standard (NM_007669)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007669
Synonyms:	CAP; CAP20; CDK; CDKI; Cdkn; Cdkn1; CI; CIP1; mda; mda6; P2; P21; p21C; p21Cip1; p21W; p21WAF; SD; SDI1; Waf; Waf1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203129) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)