

Product datasheet for **MK200163**

Mad2l1bp Mouse qPCR Template Standard (NM_025649)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_025649
Synonyms:	0610009D16Rik; 2510031P20Rik; AU018945; AW550977; C76045; CMT2A; CMT2B; Mad2lbp; mKIAA0110
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP207742) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)