

## Product datasheet for **MK200196**

### Qdpr Mouse qPCR Template Standard (NM\_024236)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_024236
Synonyms:	2610008L04Rik; D5Ertd371e; Dhpr; PKU2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212279) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)