

## Product datasheet for **MK200204**

### Rassf1 Mouse qPCR Template Standard (NM\_019713)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019713
Synonyms:	123F2; AA536941; AU044980; NORE2A; Rassf1A; Rassf1B; Rassf1C; RDA32; REH3P21
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212444) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)