

Product datasheet for **MK200221**

Arfgap3 Mouse qPCR Template Standard (NM_025445)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_025445

Synonyms: 0610009H19Rik; 1810004P07Rik; 1810035F16Rik; 9130416J18Rik; Arfgap1

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP200421) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)