

Product datasheet for **MK200222**

Gnpda2 Mouse qPCR Template Standard (NM_001038015)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001038015
Synonyms:	4921523I18Rik; 4933412A11Rik; BB189630; Gnp2; GNPDA; mKIAA4008; Sb52
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP205490) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)