

Product datasheet for **MK200224**

Gzmb Mouse qPCR Template Standard (NM_013542)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_013542
Synonyms:	AI553453; CCP; CCP-1/C1; CCP-1/C11; CCP1; Ctl; Ctla; Ctla-1; Ctla1; GZB
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP205737) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)