

Product datasheet for **MK200239**

Sike1 Mouse qPCR Template Standard (NM_025679)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_025679

Synonyms: 2810005O12Rik; 5730470L24Rik; AI450236; AI839862; Sike; Sikeb

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP214091) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)