

Product datasheet for **MK200249**

Anapc15 Mouse qPCR Template Standard (NM_027532)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_027532

Synonyms: 3200002M19Rik; 6330414C15Rik; APC15

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP213700) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)