

Product datasheet for **MK200283**

Slc40a1 Mouse qPCR Template Standard (NM_016917)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_016917
Synonyms:	Dusg; Fpn1; IREG1; MTP; MTP1; OI5; Pcm; Slc11a3; Slc39a1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP215884) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)