

## Product datasheet for **MK200290**

### **Atp6v1e1 Mouse qPCR Template Standard (NM\_007510)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_007510
<b>Synonyms:</b>	2410029D23Rik; Atp6e; Atp6e2; Atp6v1e; D6ErtD385e; E2; P31; Vma4
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP218192) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)