

Product datasheet for **MK200317**

Clip1 Mouse qPCR Template Standard (NM_019765)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019765
Synonyms:	1110007112Rik; 4631429H07Rik; AV017631; C81039; CLIP-170; Clip50; CLIP170; CYLN1; mKIAA4046
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201677) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)