

Product datasheet for **MK200325**

Dab2 Mouse qPCR Template Standard (NM_001037905)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_001037905

Synonyms: 5730435J12Rik; AA960054; AI957090; D15Wsu122e; D630005B22Rik; Doc-2; Doc2; p96

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP203710) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)