

Product datasheet for **MK200375**

Ube2h Mouse qPCR Template Standard (NM_009459)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_009459

Synonyms: 1500009C23Rik; AI181839; AI326965; AI462483; AW227540; E2-20K; Gid3; N28148; Ubc8

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP217946) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C .



[View online »](#)