

Product datasheet for **MK200379**

Spi1 Mouse qPCR Template Standard (NM_011355)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_011355
Synonyms:	Dis-1; Dis1; PU.1; Sfpi-1; Sfpi1; Spi-1; Tcfpu1; Tfpu.1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP215337) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)