

## Product datasheet for **MK200419**

### **Kcnn4 Mouse qPCR Template Standard (NM\_008433)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_008433
<b>Synonyms:</b>	IK1; IKCA1; KCa3.1; KCA4; mIKCa1; SK4; SKCas
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP211285) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)