

Product datasheet for **MK200435**

Gabpb1 Mouse qPCR Template Standard (NM_207669)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_207669
Synonyms:	BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4Tf1B; GABPB; GABPB-1; GABPB-2; GABPB1-1; GABPB1-2; NRF2B1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP205178) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)