

Product datasheet for **MK200443**

Atp6v1h Mouse qPCR Template Standard (NM_133826)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_133826
Synonyms:	0710001F19Rik; AU022349; CGI-11; SFD; SFDalpha; SFDbeta; VMA13
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201105) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)