

## Product datasheet for **MK200445**

### H2ax Mouse qPCR Template Standard (NM\_010436)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_010436

**Synonyms:** AW228881; gammaH; gammaH2ax; H2A.; H2A.X; H2af; H2afx; Hist5-; Hist5-2ax

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP205900) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)