

Product datasheet for **MK200453**

Polr3f Mouse qPCR Template Standard (NM_029763)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_029763
Synonyms:	2810411G20Rik; 3010019O03Rik; 3110032A07Rik; RPC6; RPC39
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP211223) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)