

Product datasheet for **MK200460**

B4galt3 Mouse qPCR Template Standard (NM_020579)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_020579
Synonyms:	9530061M23Rik; AA104562; AW125175; ESTM6; ESTM26; R74981
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP218014) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)