

Product datasheet for **MK200500**

Lims1 Mouse qPCR Template Standard (NM_026148)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026148
Synonyms:	2310016J22Rik; 4921524A02Rik; AI507642; AU021743; AW551584; C430041B13Rik; Lims1l; PINCH; PINCH1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP207511) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)