

## Product datasheet for **MK200670**

### **Ddx56 Mouse qPCR Template Standard (NM\_026538)**

#### **Product data:**

|                         |                                                                                                                                                                                                                                                                                               |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>    | qPCR Template Standards                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>            | NM_026538                                                                                                                                                                                                                                                                                     |
| <b>Synonyms:</b>        | 2600001H07Rik; D11Ert619e; Noh61                                                                                                                                                                                                                                                              |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (MP203413) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT                                                                                                                                                                                                                      |
| <b>Storage:</b>         | Store at -20°C.                                                                                                                                                                                                                                                                               |
| <b>Stability:</b>       | The primer mix and template standard are stable for one year from date of shipping.                                                                                                                                                                                                           |

