

## Product datasheet for **MK200688**

### Ogdh Mouse qPCR Template Standard (NM\_010956)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010956
Synonyms:	2210403E04Rik; 2210412K19Rik; AA409584; d1401; mKIAA4192
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP210504) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)