

Product datasheet for **MK200689**

Hspb8 Mouse qPCR Template Standard (NM_030704)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_030704
Synonyms:	AU018630; AW413033; Cryac; D5Ucla4; E2IG1; H11; H11K; HSP20-like; HSP22
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP205963) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)