

Product datasheet for **MK200719**

Pdzk1 Mouse qPCR Template Standard (NM_021517)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_021517

Synonyms: 1700023D20Rik; 2610507N21Rik; 4921513F16Rik; AI267131; AI314638; AL022680; D3ErtD537e

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP210652) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C .



[View online »](#)