

Product datasheet for **MK200903**

H2-Aa Mouse qPCR Template Standard (NM_010378)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010378
Synonyms:	Aalpha; H-2Aa; H2Aa; I-Aalpha; Ia-1; Ia1; IAalpha
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP206133) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)