

## Product datasheet for **MK201004**

### **Rps21 Mouse qPCR Template Standard (NM\_025587)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**ACCN:** NM\_025587

**Synonyms:** 1810049N11Rik; 2410030A14Rik

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP212915) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard:  $50 \times 10^7$  copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)