

Product datasheet for **MK201048**

Myl12b Mouse qPCR Template Standard (NM_023402)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_023402
Synonyms:	1500001M02Rik; C77744; Mylc2b; RLC-B
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP208548) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)