

Product datasheet for **MK201073**

Eif2s1 Mouse qPCR Template Standard (NM_026114)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_026114

Synonyms: 35kDa; 0910001O23Rik; 2410026C18Rik; Eif2a; eIF2alpha

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP204384) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)