

## Product datasheet for **MK201083**

### **Paccin1 Mouse qPCR Template Standard (NM\_178365)**

#### **Product data:**

|                  |                                                                                                                                                                                                                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | qPCR Template Standards                                                                                                                                                                                                                                                                       |
| ACCN:            | NM_178365                                                                                                                                                                                                                                                                                     |
| Synonyms:        | A830061D09Rik; H74; mKIAA1379; syndapin                                                                                                                                                                                                                                                       |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (MP211927) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT                                                                                                                                                                                                                      |
| Storage:         | Store at -20°C.                                                                                                                                                                                                                                                                               |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.                                                                                                                                                                                                           |



[View online »](#)