

## Product datasheet for **MK201119**

### **Nme4 Mouse qPCR Template Standard (NM\_019731)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**ACCN:** NM\_019731

**Synonyms:** 2610027N22Rik; 2810024O08Rik; 5730493H09Rik; NM23-M4; Nm23M4

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP209009) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)