

Product datasheet for **MK201138**

Pdcd6ip Mouse qPCR Template Standard (NM_011052)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_011052
Synonyms:	AI480591; Aip1; Alix; AW544830; C76364; Eig2; mKIAA1375
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP211712) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)