

Product datasheet for **MK201298**

Apoa1 Mouse qPCR Template Standard (NM_009692)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_009692
Synonyms:	Al; Alp-1; Ap; apo-Al; Apoa-1; apoA-I; Brp-; Brp-14; Ltw-; Ltw-1; Lvtw; Lvtw-1; Se; Sep; Sep-1; Sep-2; Sep2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP200800) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)