

## Product datasheet for **MK201299**

### Top3b Mouse qPCR Template Standard (NM\_011624)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | NM_011624   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (MP217040) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at $-20^{\circ}\text{C}$ .  |



[View online »](#)