

Product datasheet for **MK201319**

Usp45 Mouse qPCR Template Standard (NM_152825)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_152825
Synonyms:	3110003C05Rik; 4930550B20Rik; AI843191; Gcap7
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217921) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)