

## Product datasheet for **MK201367**

### Hopx Mouse qPCR Template Standard (NM\_175606)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_175606
Synonyms:	1110018K11Rik; 1200015P04Rik; 2300002F06Rik; AI848177; AW490897; Cameo; Hdop; Hod; Hop
Components:	1 vial of lyophilized SybGREEN qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)