

## Product datasheet for **MK201373**

### Cuta Mouse qPCR Template Standard (NM\_026307)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026307
Synonyms:	0610039D01Rik; 1810022E02Rik; 1810060C03Rik; 2700094G22Rik; AI326454
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203069) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)