

Product datasheet for **MK201379**

Fkbp10 Mouse qPCR Template Standard (NM_010221)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010221
Synonyms:	AI325255; FKBP-10; FKBP-65; Fkbp-rs1; Fkbp1-rs; Fkbp6; FKBP65; Fkbprp
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP204870) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)