

Product datasheet for **MK201396**

Fcer1g Mouse qPCR Template Standard (NM_010185)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010185
Synonyms:	AI573376; CD23; Fce1g; FcepsilonRI; FcR-gamma; FcRgamma; FcR[g]; Ly-50
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP204734) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)