

## Product datasheet for **MK201397**

### **Fcer1g Mouse qPCR Template Standard (NM\_010185)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**ACCN:** NM\_010185

**Synonyms:** AI573376; CD23; Fce1g; FcepsilonRI; FcR-gamma; FcRgamma; FcR[g]; Ly-50

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP204734) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)