

## Product datasheet for **MK201421**

### PPP2R2B Mouse qPCR Template Standard (NM\_028392)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_028392

**Synonyms:** 2900026H06Rik; 6330404L05Rik; E130009M08Rik; PP2A-PR55B; PR55-BETA; SCA12

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP220748) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard:  $50 \times 10^7$  copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at  $-20^{\circ}\text{C}$ .



[View online »](#)