

Product datasheet for **MK201437**

Hdc Mouse qPCR Template Standard (NM_008230)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_008230

Synonyms: AW108189; Hdc-a; Hdc-c; Hdc-e; Hdc-s

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP206122) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)