

## Product datasheet for **MK201438**

### **Nsdhl Mouse qPCR Template Standard (NM\_010941)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_010941
<b>Synonyms:</b>	AI747449; Bpa; H105E3; Str; XAP104
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP208645) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)