

Product datasheet for **MK201468**

PD-L1 (Cd274) Mouse qPCR Template Standard (NM_021893)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021893
Synonyms:	A530045L16Rik; B7h1; PD-; Pdcd1l; Pdcd1l1; Pdcd1lg1; Pdl1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201906) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)